



**Flexineb® is an easy-to-use equine nebulizer that silently produces a very fine mist of aerosolized medication for the fast, safe and efficient treatment of respiratory conditions in horses.**

**To help manage RAO, IAD,  
— EIPH & Other Respiratory  
Issues:**

**Antibiotics**

**Bronchodilators**

**Corticosteroids**

**Mucolytics**

**Natural Oils**

**Saline Solution**



Around the world, horses of all disciplines, breeds and ages are hindered by respiratory problems, allergies, airway irritation and inflammation. While some horses are genetically predisposed to these issues, others are impacted by changes in the season, frequent travel or a pesky allergy.

The Flexineb equine nebulizer offers a silent, portable and easy-to-use alternative for managing these respiratory issues through a highly-efficient nebulization process with fewer side effects than traditional treatments.

Flexineb® North America

JS Equine Division

[www.flexineb.us](http://www.flexineb.us)

Toll-Free

1.800.525.4339

Phone

1.731.885.6690

Fax

1.731.885.6080

Email

[info@flexineb.us](mailto:info@flexineb.us)



# Why Flexineb?



## SIMPLE & INTEGRATED

Easy setup and operation with no hoses, cords or noisy compressors.

## RAPID

10mL of saline aerosolized in 10 minutes.

## ADAPTABLE

May be used with a wide variety of medications and with metered dose inhalers.

## SILENT

Truly silent unlike most nebulizers.

## SAVINGS

On drug administration.

In order for the aerosolized medication to reach the lower airways for the treatment of respiratory issues, particle size must be less than five microns in diameter. Thanks to Flexineb's state-of-the-art components and unparalleled design, 71% of particles nebulized by the Flexineb® reaches the lower airway.

## WHY PARTICLE SIZE MATTERS



In a study conducted in conjunction with Hagyard Equine Medical Institute under the direction of Dr. Nathan Slovis in Lexington, Kentucky, the Flexineb unit nebulized radiopharmaceuticals, which were inhaled directly into the lungs. Using nuclear imagery, the dispersal of these particles via nebulization can clearly be seen to penetrate deep into the lower airway system.



Flexineb® North America

JS Equine Division

[www.flexineb.us](http://www.flexineb.us)



Nortev.com | Product of Ireland

Toll-Free

1.800.525.4339

Phone

1.731.885.6690

Fax

1.731.885.6080

Email

[info@flexineb.us](mailto:info@flexineb.us)